


Product Datasheet

ACO2 antibody (Biotin) (orb389353)

Description	Goat polyclonal antibody against ACO2 conjugated to Biotin.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Target	Aconitase 2, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 90kDa band observed in Human Heart lysates (calculated MW of 85.4kDa according to NP_001089.1).
Protein Sequence	QHVDVSPTSQRLQ.
Clonality	Polyclonal
MW	85.4
Entrez	50
Dilution Range	ELISA: 1:2000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...