

Product Datasheet

PACSIN1 antibody (orb38940)

Description Purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC-P, WB

Immunogen Synthetic Peptide

Target PACSIN1

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 50966

Uniprot ID [Q9BY11](#)


NCBI [NP_001186512.1](#), [NP_065855.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of HL-60 cell line...