


## Product Datasheet

### Mouse KLK1/Kallikrein 1 ELISA Kit (orb389535)

<b>Description</b>	Mouse KLK1/Kallikrein 1 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	62.5 pg/ml - 4,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: I25-D261
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	Dilution ratio of 1:100, concentration in serum and plasma is 50-80 ng/ml
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
<b>MW</b>	28775 MW
<b>Uniprot ID</b>	<b>P15947</b>
<b>Sensitivity</b>	12 pg/ml
<b>Expiration Date</b>	Please enquire.



Mouse KLK1  
PicoKine  
ELISA Kit  
standard c...