

Product Datasheet

RPN1 antibody (orb390123)

Description Rabbit polyclonal antibody to RPN1.

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Recombinant fusion protein of human RPN1(NP_002941.1).

Preservatives PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 68kDa


Uniprot ID **P04843**

Dilution Range WB: 1:500-2000, IHC-P: 1:50-200


Expiration Date 12 months from date of receipt.



Immunohistochemical staining of human colon tissue.



Immunohistochemical staining of human liver tissue.



Immunohistochemical staining of human tissue.