



Product Datasheet

GABARAP antibody (orb390149)


| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to GABARAP. |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC, WB |
| Immunogen | A synthetic peptide of human GABARAP (NP_009209.1). |
| Preservatives | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Affinity purification |
| MW | 14kDa |
| Uniprot ID | O95166 |
| Entrez | 11337 |
| Dilution Range | WB: 1:500-2000, IHC-P: 1:50-200 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical staining of mouse ki...



Immunohistochemical staining of human li...



Immunohistochemical staining of rat hear...