

## Product Datasheet

### GM13125 antibody (orb390170)

**Description** Rabbit polyclonal antibody to GM13125.

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** IHC, WB

**Immunogen** Recombinant fusion protein of mouse Gm13125(NP\_001108549.1).

**Preservatives** PBS with 0.02% sodium azide,50% glycerol,pH7.3.

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Isotype** IgG


**Clonality** Polyclonal

**MW** 55kDa


**Uniprot ID** [Q5VT98](#)

**Dilution Range** WB: 1:500-2000, IHC-P: 1:50-200


**Expiration Date** 12 months from date of receipt.



Immunohistochemical staining of human st...



Immunohistochemical staining of human pr...



Immunohistochemical staining of human li...