



Product Datasheet

SCARB2 antibody (orb390190)


Description	Rabbit polyclonal antibody to SCARB2.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-380 of human SR-B2/LIMP2 (NP_005497.1).
Preservatives	PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80kDa
Uniprot ID	Q14108
Dilution Range	WB: 1:1000-3000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of L929 cell...



Immunofluorescence analysis of U2OS cell...



Western blot analysis of various lysates...