

## Product Datasheet

### EIF3J antibody (orb390316)

**Description** Rabbit polyclonal antibody to EIF3J.

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** ICC, IF, WB

**Immunogen** Recombinant fusion protein of human EIF3J(NP\_003749.2).

**Preservatives** PBS with 0.01% thiomersal,50% glycerol,pH7.3.

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note** For research use only

**Isotype** IgG


**Clonality** Polyclonal

**MW** 37kDa


**Uniprot ID** [O75822](#)

**Dilution Range** WB: 1:200-2000


**Expiration Date** 12 months from date of receipt.



Immunohistochemical staining of mouse sp...



Immunohistochemical staining of human ga...



Immunohistochemical staining of human lu...