


Product Datasheet

BRCA1 (Phospho-S988) antibody (orb390335)

Description	Rabbit polyclonal antibody to BRCA1 (Phospho-S988).
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	A synthetic phosphorylated peptide around S988 of human BRCA1(NP_009225.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	7kDa/78-85kDa/202-210kDa
Uniprot ID	P38398
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...