

## Product Datasheet

# Phenylalanine ammonia-lyase Microplate Assay Kit (orb390724)

<b>Description</b>	Detection and Quantification of Phenylalanine ammonia-lyase Activity.
<b>Reactivity</b>	All
<b>Tested Applications</b>	FA
<b>Storage</b>	Shipped and store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
<b>Sample Types</b>	Tissue extracts, Cell lysate
<b>CAS Number</b>	9024-28-6
<b>Expiration Date</b>	12 months from date of receipt.