



Product Datasheet

Glutathione Reductase Microplate Assay Kit (orb390749)

Description	Detection and Quantification of Glutathione Reductase Activity.
Reactivity	All
Range	4 µmol/L - 400 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-48-3
Expiration Date	12 months from date of receipt.



Representative Standard Curve



Positive Control reaction in 96-well pla...