

Product Datasheet

NADH Oxidase Microplate Assay Kit (orb390750)

Description	Detection and Quantification of NADH Oxidase Activity.
Reactivity	All
Range	30 $\mu\text{mol/L}$ - 600 $\mu\text{mol/L}$
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
CAS Number	9032-21-7
Expiration Date	12 months from date of receipt.