


Product Datasheet

Ascorbate Peroxidase Microplate Assay Kit (orb390758)

Description	Detection and Quantification of Ascorbate Peroxidase Activity.
Reactivity	All
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
Expiration Date	12 months from date of receipt.



Samples were assayed using the 96-well m...