

Product Datasheet

Ceruloplasmin Microplate Assay Kit (orb390765)

Description	Detection and Quantification of Ceruloplasmin Activity.
Reactivity	All
Range	0.01 U - 1 U
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Serum, Plasma, Other biological fluids
Sample Types	Serum, Plasma, Other biological fluids
CAS Number	9031-37-2
Expiration Date	12 months from date of receipt.