


Product Datasheet

Tannin Microplate Assay Kit (orb390766)

Description	Detection and Quantification of Tannin Content.
Reactivity	All
Range	10 µg/ml - 1000 µg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
CAS Number	1401-55-4
Expiration Date	12 months from date of receipt.



Representative Standard Curve