

## Product Datasheet

### Peroxidase Microplate Assay Kit (orb390768)

**Description** Detection and Quantification of Peroxidase Activity.

**Reactivity** All

**Tested Applications** FA


**Storage** Shipped and store at 4°C for 6 months.

**Note** For research use only

**Application notes** Detection Method: Colorimetric  
Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids

**Sample Types** Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

**Expiration Date** 12 months from date of receipt.



Positive Control reaction in 96-well pla...