


Product Datasheet

Glutaminase Microplate Assay Kit (orb390771)

Description	Detection and Quantification of Glutaminase Activity.
Reactivity	All
Range	0.01 mmol/L - 1 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate
Sample Types	Tissue extracts, Cell lysate
CAS Number	9001-47-2
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve