

Product Datasheet

Pyruvate Kinase Microplate Assay Kit (orb390780)

Description Detection and Quantification of Pyruvate Kinase Activity.

Reactivity All

Range 4 µmol/L - 400 µmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids


Sample Types Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 9001-59-6

Expiration Date 12 months from date of receipt.



Representative Standard Curve



Positive Control reaction in 96-well pla...