

## Product Datasheet

# Phosphoenolpyruvate Carboxylase Microplate Assay Kit (orb390782)

**Description** Detection and Quantification of Phosphoenolpyruvate Carboxylase Activity.

**Reactivity** All

**Range** 4 µmol/L - 400 µmol/L

**Tested Applications** FA

**Storage** Shipped and store at 4°C for 6 months.


**Note** For research use only

**Application notes** Detection Method: Colorimetric  
Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids


**Sample Types** Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

**CAS Number** 9067-77-0

**Expiration Date** 12 months from date of receipt.



Representative Standard Curve



Positive Control reaction in 96-well pla...