

Product Datasheet

Neutral Protease Microplate Assay Kit (orb390784)

Description Detection and Quantification of Neutral Protease Activity.

Reactivity All

Range 0.05 mmol/L - 5 mmol/L

Tested Applications FA

Storage Shipped and store at 4°C for 6 months.


Note For research use only

Application notes Detection Method: Colorimetric
Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids


Sample Types Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

CAS Number 9068-59-1

Expiration Date 6 months at 4°C



Representative Standard Curve



Positive Control reaction in 96-well pla...