



## Product Datasheet

# Pyruvate Decarboxylase Microplate Assay Kit (orb390791)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Detection and Quantification of Alcohol Dehydrogenase Activity.  |
| <b>Reactivity</b>          | All  |
| <b>Range</b>               | 4 µmol/L - 400 µmol/L  |
| <b>Tested Applications</b> | FA   |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months.   |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detection Method: Colorimetric<br>Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids |
| <b>Sample Types</b>        | Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids  |
| <b>CAS Number</b>          | 9001-04-1  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Representative Standard Curve



Positive Control reaction in 96-well pla...