


## Product Datasheet

### Pectinase Microplate Assay Kit (orb390793)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Detection and Quantification of Pectinase Activity.  |
| <b>Reactivity</b>          | All  |
| <b>Range</b>               | 0.1 mg/ml - 4 mg/ml  |
| <b>Tested Applications</b> | FA   |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months.   |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detection Method: Colorimetric<br>Sample Types: Tissue extracts, Cell lysate, Cell culture media and Other biological fluids |
| <b>Sample Types</b>        | Tissue extracts, Cell lysate, Cell culture media, Other biological fluids  |
| <b>CAS Number</b>          | 9032-75-1  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |



Representative Standard Curve

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)