


Product Datasheet

Pectinase Microplate Assay Kit (orb390793)

Description	Detection and Quantification of Pectinase Activity.
Reactivity	All
Range	0.1 mg/ml - 4 mg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9032-75-1
Expiration Date	12 months from date of receipt.



Representative Standard Curve