


Product Datasheet

Chitinase Microplate Assay Kit (orb390794)

Description	Detection and Quantification of Chitinase Activity.
Reactivity	All
Range	0.01 mg/ml - 1 mg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-06-3
Expiration Date	12 months from date of receipt.



Representative Standard Curve