



Product Datasheet

Neutral Xylanase Microplate Assay Kit (orb390795)

Description	Detection and Quantification of Neutral Xylanase Activity.
Reactivity	All
Range	0.01 mmol/L - 1 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Animal feeds, Enzyme preparations, Bread improver mixtures and other materials
Sample Types	Animal feeds, Enzyme preparations, Bread improver mixtures, other materials
CAS Number	9025-57-4
Expiration Date	12 months from date of receipt.



Representative Standard Curve



Positive Control reaction in 96-well pla...