



## Product Datasheet

### Acidic Xylanase Microplate Assay Kit (orb390796)

<b>Description</b>	Detection and Quantification of Acidic Xylanase Activity.
<b>Reactivity</b>	All
<b>Range</b>	0.01 mmol/L - 1 mmol/L
<b>Tested Applications</b>	FA
<b>Storage</b>	Shipped and store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric Sample Types: Animal feeds, Enzyme preparations, Bread improver mixtures and other materials
<b>Sample Types</b>	Animal feeds, Enzyme preparations, Bread improver mixtures, other materials
<b>CAS Number</b>	9025-57-4
<b>Expiration Date</b>	12 months from date of receipt.



Representative Standard Curve



Positive Control reaction in 96-well pla...