


Product Datasheet

Beta-xylosidase Microplate Assay Kit (orb390798)

Description	Detection and Quantification of Beta-xylosidase Activity.
Reactivity	All
Range	10 µmol/L - 1000 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate and Other biological fluids
Sample Types	Tissue extracts, Cell lysate, Other biological fluids
Expiration Date	12 months from date of receipt.



Representative Standard Curve