

Product Datasheet

Pectate Lyase Microplate Assay Kit (orb390799)

| | |
|----------------------------|--|
| Description | Detection and Quantification of Pectate Lyase Activity. |
| Reactivity | All |
| Tested Applications | FA |
| Storage | Shipped and store at 4°C for 6 months. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate and Other biological fluids |
| Sample Types | Tissue extracts, Cell lysate, Other biological fluids |
| CAS Number | 9015-75-2 |
| Expiration Date | 12 months from date of receipt. |