


Product Datasheet

Cellulase Microplate Assay Kit (orb390801)

Description	Detection and Quantification of Cellulase Activity.
Reactivity	All
Range	0.3 mmol/L - 3 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9012-54-8
Expiration Date	12 months from date of receipt.



Representative Standard Curve