

## Product Datasheet

# Beta-Glucosidase Microplate Assay Kit (orb390803)

**Description** Detection and Quantification of Beta-Glucosidase Activity.

**Reactivity** All

**Range** 10 µmol/L - 1000 µmol/L

**Tested Applications** FA

**Storage** Shipped and store at 4°C for 6 months.


**Note** For research use only

**Application notes** Detection Method: Colorimetric  
Sample Types: Tissue extracts, Cell lysate, Cell culture media and Other biological fluids


**Sample Types** Tissue extracts, Cell lysate, Cell culture media, Other biological fluids

**CAS Number** 9001-22-3

**Expiration Date** 12 months from date of receipt.



Representative Standard Curve



Positive Control reaction in 96-well pla...