



Product Datasheet

Beta-Galactosidase Microplate Assay Kit (orb390805)

Description	Detection and Quantification of Beta-Galactosidase Activity.
Reactivity	All
Range	10 µmol/L - 1000 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9031-11-2
Expiration Date	12 months from date of receipt.



Representative Standard Curve



Positive Control reaction in 96-well pla...