


Product Datasheet

Sorbitol Microplate Assay Kit (orb390806)

Description	Detection and Quantification of Sorbitol Content.
Reactivity	All
Range	0.1 mg/ml - 5 mg/ml
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	50-70-4
Expiration Date	12 months from date of receipt.



Representative Standard Curve