


Product Datasheet

Sorbitol Dehydrogenase Microplate Assay Kit (orb390807)

Description	Detection and Quantification of Sorbitol Dehydrogenase Activity.
Reactivity	All
Range	30 µmol/L - 600 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media and Other biological fluids
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9028-21-1
Expiration Date	12 months from date of receipt.



Representative Standard Curve