

## Product Datasheet

### Serum Kalium Microplate Assay Kit (orb390810)

**Description** Detection and Quantification of Serum Kalium Content.

**Reactivity** All

**Range** 0.2 mmol/L - 2 mmol/L

**Tested Applications** FA


**Storage** Shipped and store at 4°C for 6 months.

**Note** For research use only

**Application notes** Detection Method: Colorimetric  
Sample Types: Serum

**Sample Types** Serum, Urine, Saliva, Other biological fluids

**Expiration Date** 12 months from date of receipt.



Representative  
Standard  
Curve