


Product Datasheet

Serum Iron Microplate Assay Kit (orb390812)

Description	Detection and Quantification of Serum Ferrum Content.
Reactivity	All
Range	5 µmol/L - 500 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Serum
Sample Types	Serum, Plasma, Urine, Saliva, Tissue extracts, Cell lysate, Other biological fluids
Expiration Date	12 months from date of receipt.



Representative Standard Curve