


## Product Datasheet

# Serum Phosphorus Microplate Assay Kit (orb390814)

<b>Description</b>	Detection and Quantification of Serum Phosphorus Content.
<b>Reactivity</b>	All
<b>Range</b>	0.01 mmol/L - 0.4 mmol/L
<b>Tested Applications</b>	FA
<b>Storage</b>	Shipped and store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric Sample Types: Serum
<b>Sample Types</b>	Serum, Urine, Saliva, Other biological fluids
<b>Expiration Date</b>	12 months from date of receipt.



Representative  
Standard  
Curve