


Product Datasheet

Human ADIPOQ ELISA Kit (orb390822)

Description	Sandwich ELISA for quantitative detection of human ADIPOQ concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	ADIPOQ
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q15848
Sensitivity	15 pg/ml
Entrez	9370
Expiration Date	Please enquire.



Representative
Standard
Curve