


Product Datasheet

Human ANGPT1 ELISA Kit (orb390824)

| | |
|----------------------------|---|
| Description | Sandwich ELISA for quantitative detection of human ANGPT1 concentrations. |
| Reactivity | Human |
| Range | 31.2 pg/ml - 2000 pg/ml |
| Tested Applications | ELISA |
| Target | ANGPT1 |
| Storage | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma |
| Sample Types | Cell Culture Supernates, Serum, Plasma |
| Uniprot ID | Q15389 |
| Sensitivity | 4 pg/ml |
| Entrez | 284 |
| Expiration Date | Please enquire. |



Representative Standard Curve