


## Product Datasheet

### Human BMP4 ELISA Kit (orb390834)

<b>Description</b>	Sandwich ELISA for quantitative detection of human BMP4 concentrations.
<b>Reactivity</b>	Human
<b>Range</b>	62.5 pg/ml - 4000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	BMP4
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric Sample Types: Bone Tissue Homogenates, Cell Culture Supernatants
<b>Sample Types</b>	Bone Tissue Homogenates, Cell Culture Supernates
<b>Uniprot ID</b>	<b>P12644</b>
<b>Sensitivity</b>	4 pg/ml
<b>Entrez</b>	<b>652</b>
<b>Expiration Date</b>	Please enquire.



Representative  
Standard  
Curve