


Product Datasheet

Human CTSB ELISA Kit (orb390837)

Description	Sandwich ELISA for quantitative detection of human CTSB concentrations.
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	CTSB
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P07858
Sensitivity	5 pg/ml
Entrez	1508
Expiration Date	Please enquire.



Representative
Standard
Curve