


## Product Datasheet

### Human RANTES ELISA Kit (orb390852)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Sandwich ELISA for quantitative detection of human RANTES concentrations.                                    |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 7.8 pg/ml - 500 pg/ml  |
| <b>Tested Applications</b> | ELISA  |
| <b>Target</b>              | CCL5   |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.                |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detection Method: Colorimetric<br>Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma, Tissue Homogenates   |
| <b>Uniprot ID</b>          | <b>P13501</b>  |
| <b>Sensitivity</b>         | 4 pg/ml  |
| <b>Entrez</b>              | <b>6352</b>  |
| <b>Expiration Date</b>     | Please enquire.  |



Representative  
Standard  
Curve