


Product Datasheet

Human CXCL6 ELISA Kit (orb390884)

| | |
|----------------------------|---|
| Description | Sandwich ELISA for quantitative detection of human CXCL6 concentrations. |
| Reactivity | Human |
| Range | 31.2 pg/ml - 2000 pg/ml |
| Tested Applications | ELISA |
| Target | CXCL6 |
| Storage | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma |
| Sample Types | Cell Culture Supernates, Serum, Plasma |
| Uniprot ID | P80162 |
| Sensitivity | 10 pg/ml |
| Entrez | 6372 |
| Expiration Date | Please enquire. |



Representative
Standard
Curve