


Product Datasheet

Human Cystatin C ELISA Kit (orb390889)

Description	Sandwich ELISA for quantitative detection of human Cystatin C concentrations.
Reactivity	Human
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	CST3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P01034
Sensitivity	10 pg/ml
Entrez	1471
Expiration Date	Please enquire.



Representative
Standard
Curve