


## Product Datasheet

### Human CDH1 ELISA Kit (orb390893)

<b>Description</b>	Sandwich ELISA for quantitative detection of human CDH1 concentrations.
<b>Reactivity</b>	Human
<b>Range</b>	156 pg/ml - 10000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	CDH1
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine
<b>Sample Types</b>	Cell Culture Supernates, Serum, Plasma, Urine
<b>Uniprot ID</b>	<b>P12830</b>
<b>Sensitivity</b>	15 pg/ml
<b>Entrez</b>	<b>999</b>
<b>Expiration Date</b>	Please enquire.



Representative  
Standard  
Curve