


Product Datasheet

Human IL18 ELISA Kit (orb390934)

Description	Sandwich ELISA for quantitative detection of human IL18 concentrations.
Reactivity	Human
Range	39 pg/ml - 2500 pg/ml
Tested Applications	ELISA
Target	IL18
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q14116
Sensitivity	20 pg/ml
Entrez	3606
Expiration Date	Please enquire.



Representative
Standard
Curve