


Product Datasheet

Human IL18BP ELISA Kit (orb390936)

Description	Sandwich ELISA for quantitative detection of human IL18BP concentrations.
Reactivity	Human
Range	93.7 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	IL18BP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	095998
Sensitivity	10 pg/ml
Entrez	10068
Expiration Date	Please enquire.



Representative Standard Curve