


## Product Datasheet

### Human KLK6 ELISA Kit (orb390946)

<b>Description</b>	Sandwich ELISA for quantitative detection of human KLK6 concentrations.
<b>Reactivity</b>	Human
<b>Range</b>	78 pg/ml - 5000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	KLK6
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Cell Lysates, Serum, Plasma, Tissue Homogenates
<b>Sample Types</b>	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
<b>Uniprot ID</b>	<b>Q92876</b>
<b>Sensitivity</b>	3 pg/ml
<b>Entrez</b>	<b>5653</b>
<b>Expiration Date</b>	Please enquire.



Representative Standard Curve