


## Product Datasheet

### Human MIF ELISA Kit (orb390963)

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Sandwich ELISA for quantitative detection of human MIF concentrations.                                       |
| <b>Reactivity</b>          | Human                                                                                                        |
| <b>Range</b>               | 31.2 pg/ml - 2000 pg/ml                                                                                      |
| <b>Tested Applications</b> | ELISA                                                                                                        |
| <b>Target</b>              | MIF                                                                                                          |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.                |
| <b>Note</b>                | For research use only                                                                                        |
| <b>Application notes</b>   | Detection Method: Colorimetric<br>Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma, Tissue Homogenates                                                   |
| <b>Uniprot ID</b>          | <b>P14174</b>                                                                                                |
| <b>Sensitivity</b>         | 15 pg/ml                                                                                                     |
| <b>Entrez</b>              | <b>4282</b>                                                                                                  |
| <b>Expiration Date</b>     | Please enquire.                                                                                              |



Representative  
Standard  
Curve