


Product Datasheet

Human RETN ELISA Kit (orb390998)

Description	Sandwich ELISA for quantitative detection of human RETN concentrations.
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	RETN
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q9HD89
Sensitivity	4 pg/ml
Entrez	56729
Expiration Date	Please enquire.



Representative Standard Curve