


Product Datasheet

Human TGF beta 1 ELISA Kit (orb391007)

Description	Sandwich ELISA for quantitative detection of human TGF beta 1 concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	TGFB1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Urine, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates, Urine
Uniprot ID	P01137
Sensitivity	30 pg/ml
Entrez	7040
Expiration Date	Please enquire.



Representative
Standard
Curve